


# LEITZ Universal Rotating Stages and Accessories



11534-55

Universal rotating stage UT 4



9519-55

Universal rotating stage UT 5 with centring base



9508-55

Schmidt parallel slide

The universal rotating stages are supplied in two designs which differ in the number of their axes or directions of rotation and have received their designations UT 4 and UT 5. The number 4 or 5 signifies\* the number of axes. In conjunction with supplementary equipment, our rotating stages may be employed for index measurements of grains in accordance with various methods, e.g. the  $\lambda$  or  $\lambda_z$  variation, or, with the aid of the Waldmann hollow glass sphere, they may serve to carry out morphological and crystal-optical measurements on crystals.

For the optical investigation of crystals after the Fedorow method, the UT 4 and UT 5 stages fulfil all the requirements. The latter surpasses the former in that it allows of setting the second plane of symmetry immediately after the first one has been found, thus facilitating the measuring process.

All universal rotating stages have a practical device for exchanging and mounting the thin sections and allow of a simple and rapid adjustment of the section surface into the intersection of the axes of rotation. For the analysis of specimen structures it is advantageous to equip the stages UT 4 and UT 5 with the auxiliary angular slide (FEGFU) for parallel movement of the specimen which, however, necessitates a special mount of the upper segment (FEGSE).

For the demonstration of the principle underlying the construction and use of the universal rotating stages, indicatrix models of optically uniaxial and biaxial crystals can be accommodated on all types of UT-stages where they take the place of the inner stage plate.

The application of the universal rotating stage methods necessitates the use of special UM objectives which are corrected for a uniform working distance (1.5 mm.) in relation to the UT stage segment and supplied with built-in iris diaphragm. The objectives of higher numerical aperture, UM 20/0.33 and UM 32/0.30 call for a special condenser cap and are ideal for the determination of directions of reference in crystallographic work.

## Specifications

Universal rotating stage UT 5 (5 axes of rotation), on centring base, with pair of segments  $n_D = 1.554$ , in case

Universal rotating stage UT 4 (4 axes of rotation), on centring base, with pair of segments  $n_D = 1.554$ , in case

Interchangeable pair of segments  $n_D = 1.554$  (as included above)

Interchangeable pair of segments  $n_D = 1.649$

Interchangeable pair of segments  $n_D = 1.516$

Segments with other refractive index on request.

For the analysis of structures of materials with the aid of the UT 4 or UT 5 stage the following is required:

Schmidt parallel guiding slide with mm. scale, for systematically displacing the thin section under the segment

The Schmidt slide necessitates a modification to the normal mount of the upper segments, extra for each segment

For general examination of thin sections on the UT stage without upper segment:

Auxiliary stage clip

If the universal rotating stages are to be used with the microscope ORTHOLUX-POL required:

Intermediate slide with clamping screw for raising the microscope tube

ICGLI

FEDRO

KOSID

FEGNO

FEGMA

FEGFU

FEGSE

FEKLE

POKUN

## Objectives for universal stage methods

Built-in iris diaphragm, free working distance to top of segment 1.5 mm., magnification and aperture values apply to segment  $n_D = 1.554$

UM 5/0.10

UM 10/0.20

UM 20/0.33

UM 32/0.30

Auxiliary setting objective for adjustment of the  $A_4$  axis

Objective changing ring (required for each objective when the microscope to be used is equipped with an objective centring clutch)

Supplementary condenser for UM 20/0.33 and UM 32/0.30, taking the place of the normal swing-out condenser top

(a) for UT 5

(b) for UT 4

ICNXI

ICOZI

ICPBI

ICQDI

PEHOJ

PIZUT

IITLC

IITNG



9508-55

9511-55

9512-55

Auxiliary objective Supplementary condenser

## LEITZ Universal Stage Coscroscope

Special equipment has been designed to supplement any model of the universal rotating stage for conoscopic work. The following items are required to build up a universal stage coscroscope:



Condenser on dovetail slide with sleeve taking the polarizer (for use in place of the normal microscope condenser) and objective UMK 32/0.60  
 Upper segment with small hemisphere (11.5 mm. dia.) refractive index 1.516  
 Upper segment with small hemisphere (11.5 mm. dia.) refractive index 1.554  
 Upper segment with small hemisphere (11.5 mm. dia.) refractive index 1.649  
 The lower segment is the same as that supplied for the UT stages:  
 Lower segment with hemisphere 25 mm. dia., refractive index 1.516  
 Lower segment with hemisphere 25 mm. dia., refractive index 1.554  
 Lower segment with hemisphere 25 mm. dia., refractive index 1.649  
 Objective UMK 32/0.60  
 Objective UMK 50/0.60  
 (Special system for the conoscopic examination of small objects)

KOSOF  
 KOREB  
 KORIC  
 KOROD

KORUF  
 KOSAP  
 KOSEC  
 IIWUS  
 PEHIX

## The Waldmann Hollow Glass Sphere



The Waldmann hollow glass sphere

The Waldmann hollow glass sphere in position on the UT 4 stage. Case with all the accessories at the right



This accessory to the universal rotating stages which can be mounted on our polarizing microscopes is for the morphological and crystal-optical examination of crystals ranging in diameter from 1 to 11 mm. It consists of a hollow sphere of optical glass, 27 mm. in diameter and with a 12 mm. bore. The space inside this sphere is filled up free from air bubbles with a suitable immersion fluid. The closure cap, with the object holder in place, seals the sphere without extending beyond its surface at any point. This renders the sphere capable of being turned under the microscope without limitation in every direction.

Advantages of the hollow glass sphere:

- Transparent crystals up to the stated maximum size can easily be brought into the centre of the sphere, where they can be examined without the risk of damage, and without preparatory measures.
- In the examination of thin sections, as has hitherto been customary, an initial position unfavourable for many components of the section is given by the plane of the section. By contrast, the sphere with its unlimited range of rotation allows a favourable initial position of the grain to be chosen, in addition to which the object can also be transposed on its holder.
- The angular space remaining accessible for observation in a plane of symmetry of the cap (about  $26^\circ$ ), is considerably larger than in ordinary segments, in which a spherical belt of less than  $90^\circ$  only can be fully utilized for transmitted light microscopy.
- The conoscopic examination method can also be employed at any time.

Three different types of closure caps are available for the Waldmann hollow glass sphere:




Left: closure cap with fixed glass pin  
Centre: closure cap with radially movable glass pin  
Right: closure cap with cross-pincers

1. Closure cap with fixed glass pin on which the object to be examined is cemented.
2. Closure cap with radially movable glass pin allowing the object to be brought into the centre of the hollow sphere.
3. Closure cap with cross pincers for crystals of 5–11 mm. dia. The pincers are opened and closed by means of a key.

#### Specification:

Waldmann hollow glass sphere, with tongs, adapter ring, and holder for use on UT stage, key for vertical adjustment of closure caps. 2 wooden rings as supports for the sphere, centring gauge, and 3 closure caps, in case

IRUX



Emmons circulation cell


### Emmons Circulation Cell

This supplementary equipment is primarily offered for the determination of refractive indices after the double variation method ( $\lambda, t$  method as described by Emmons).

The equipment consists of a special lower segment for UT stages and a circulation cell surrounding this segment for varying the temperature. A thermostat is also provided for the purpose. For  $\lambda$  variation when carrying out the  $\lambda$  method the linear mirror monochromator must be employed.

Lower segment with circulation cell, in case

PEJAG



Angle-true net ruling

Surface-true net ruling

### Accessories for Universal Stage Methods

For the evaluation of the measurements obtained with the aid of the universal stages the following accessories are available:

Angle-true stereographic net ruling (according to Wulff) with rotating device (according to M. Reinhard) for the tracing paper

IVWXI

Angle-true stereographic net ruling, single sheet

IZWLI

Surface-true net ruling (according to Lambert) with rotating device (according to M. Reinhard) for the tracing paper

IVXZI

Surface-true net ruling, single sheet

IZMYI


For the demonstration and explanation of the universal stage methods and the operations involved we supply, as suggested by M. Reinhard, the following indicatrix model of crystals which are accommodated on the universal stage after the stage inset has been removed:

Indicatrix model of an optically biaxial crystal

ICSHI


Indicatrix model of an optically uniaxial crystal

ICTKI



Indicatrix model of an optically biaxial crystal

Indicatrix model of an optically uniaxial crystal




9523-55

Auxiliary arrangement INDEX



9520-55

Heating ring for the hemisphere



9522-55

Colour filter



Plano-parallel cover plates



9522-55

special eyepiece INDEX 8 x

### Universal Stage Refractometer

This equipment, also known as the Berek microscope refractometer, is designed for the determination of the refractive index of grain preparations after the embedding method and with the aid of a universal stage.

The refractive index of the immersion is adapted to that of the grain by varying the temperature. After a changeover from transmitted to diffused incident light the refractive index results from the setting of total reflection on the UT-stage.

#### Specification:

Auxiliary arrangement INDEX comprising a lower segment, upper hemisphere with cavity and 3 plano-parallel cover plates, heating ring for the hemisphere with connecting cable, regulating transformer type RT 35 for 110/220 volts a. c.

Special eyepiece INDEX 8x (30 mm. dia.) with cross lines in helical focusing mount and adjustable eyelens

Objective UM 20/0.33

Illuminating stand with opal glass plate and tubular lamp 25 watts for direct mains connection\*, fitting to the microscope foot

IDCBI

IDEFI

ICPBI

PEJEH

When using a filament lamp as source of light one of the following filters must be placed on the INDEX eyepiece or the top analyser:

Orange filter

Colour filter 670

Colour filter 550

Colour filter 480

IDKQI

IDESI

IDUMI

IDNWI

\*) please specify mains voltage when ordering

Design subject to alterations without notice


**ERNST LEITZ GMBH WETZLAR GERMANY**  
 Subsidiary: Ernst Leitz (Canada) Ltd., Midland, Ontario

List 55-10b/Engl.

Printed in Germany

IV/62/DX/L










## Procedure for Examinations with the Universal Rotating Stage

### I. Preparations

1. (a) Crossing the Nicols.  
(b) Centring the objectives.
2. Estimation of the approximate degree of light refraction of the mineral to be examined as compared with Canada balsam, to permit correct selection of the appropriate segment.
3. Secure the object slide (with the rock section on top) and the spherical segments on the universal rotating stage. The object slide should be between 0.9 and 1.1 mm. thick. A little immersion oil is applied to the lower surface of the object slide and to the cover glass of the specimen; the slide is then placed on the glass plate, which rests loosely in the U-stage, in such a manner that it lies at right angles to the line connecting the two screws fastening the upper segment; these two screws press the upper segment onto the specimen with an elastic pressure. The fastening screw marked with a ring can be screwed in beforehand. The upper segment, which is designed with a corresponding slot in its mounting, is pushed under the spring-mounted socket of the screw, and is finally screwed tight with the second screw. In order to be able to remove the specimen and the upper segment conveniently, this screw is removed, and the screw marked with a ring is then loosened. A little immersion oil is also applied to the lower segment; it is then pressed against the glass plate from beneath, and is held in this position by adhesion.
4. There are 5 different axes of rotation; their varying nomenclature in the literature is listed below. Berek's nomenclature of the axes, i. e.  $A_1$ ,  $A_2$ ,  $A_3$ ,  $A_4$ ,  $A_5$  corresponds to the sequence of the required operations, from the centre outwards.  $A_1$ ,  $A_3$  and  $A_5$  (odd indices) are vertical axes of rotation, while  $A_2$  and  $A_4$  (even indices) are horizontal.

Cf. M. Berek „Neue Wege zur Universalmethode“ (New Applications of the Universal Method) N. Jb. Mineral. Vol. 48 (1923), pp. 34–62.





Berek (1924)	Nomenclature according to Nikitin-Duparc-Reinhard	(1931) Reinhard	R. C. Emmons
A <sub>1</sub>	N (normal axis)	N ( )	I. V.
A <sub>2</sub>	H (horizontal axis)	H ( )	N-S
A <sub>3</sub>	M (mobile axis)	A (auxiliary axis)	O. V.
A <sub>4</sub>	I (immobile axis)	K (control axis)	O. E-W
A <sub>5</sub>		M (Microscope axis)	M

**5. Assembly and adjustment of the U-stage on the object stage of the microscope.**

- Remove ring plate from the object stage.
- Fasten universal rotating stage to the object stage with two screws, so that the graduated drum of the A<sub>4</sub> axis is on the right.
- Lock A<sub>3</sub>; lock A<sub>2</sub> and A<sub>4</sub> at index "0".  
(The centring of the objective used can also be carried out at this point by rotating on axis A<sub>5</sub> and observing the specimen).
- Adjusting the A<sub>1</sub> or A<sub>3</sub> axis. Rotate on axis A<sub>1</sub> or A<sub>3</sub>. If the axis of rotating of the specimen does not coincide with the point of intersection of the cross-lines, then the two screws securing the U-stage are loosened, and the latter is repositioned accordingly. Centring is correct if the axis of rotation of the specimen still coincides with the point of intersection of the cross-lines on rotating on axis A<sub>5</sub>, or A<sub>1</sub> and A<sub>3</sub> respectively. The A<sub>4</sub> axis is now adjusted parallel to the horizontal line of the cross-lines of the eyepiece.
- Focus the microscope on the surface of the upper segment. Set the index of A<sub>3</sub> to 270°. Rotate on A<sub>4</sub> and observe whether the dust particles on the surface of the upper segment move exactly parallel to the vertical line of the eyepiece cross-lines; if not, rotate on A<sub>3</sub> until the required effect is achieved. This adjustment must be carried out carefully. The final position of A<sub>3</sub> is the zero position of A<sub>5</sub>. Lock, and note down this position.

- (f) Turn all axes to the normal position and lock.  
 $A_1 : 0$     $A_2 : 0$     $A_3 : 270$  (or 90)    $A_4 : 0$   
 $A_5$  : zero position as determined above.
- (g) Adjusting the vertical position of the specimen.  
 Rotate on  $A_4$  in both directions from the zero position. If on doing so, the plane of the specimen appears to move above the horizontal line of the cross-lines, then the plane of the specimen lies under the axis  $A_4$ , and must be raised. This is done with the aid of the ring fitted with four notches, which is located under the centre of the stage; this ring is turned in a clockwise direction, holding the screws securing the upper segment in position at the same time, until the axis  $A_4$  finally lies in the plane of the specimen. Conversely it is also possible to hold the notches in position, and to rotate on  $A_4$ .
- (h) Close the objective iris diaphragm. Focus the Bertrand lens on same and centre. Centre source of light to the objective iris diaphragm.
- (i) For the graphic evaluation of the results of the measurements, we supply **angle-true** (stereographic, Wulff's, Reinhard's) and **surface-true** (Lambert's, Schmidt's) net rulings with a revolving device after M. Reinhard.
- (k) It should be noted that when using the Wulff's net ruling, the index of the  $A_3$  axis lies at  $270^\circ$ , where it normally remains arrested, since the graduation of the net ruling is usually marked from the top (North) with  $0 : 360^\circ$ , and lies correspondingly at  $270^\circ$  analogous as with the U-stage.
6. In the 5-axis universal rotating stage, an additional horizontal axis has been introduced between axes  $A_2$  and  $A_3$ , and lies vertical to  $A_2$ . This axis, which bears the designation  $A'_2$  and has its angle graduation on adjustable arcs, facilitates determination of the position of the indicatrix in the crystal. The advantages of the 5-axis stage in determining the second and third planes of symmetry are described in detail in the book by Rinne-Berek: „Anleitung zu optischen Untersuchungen mit dem Polarisationsmikroskop“ (Instructions for Optical Examinations with the Polarizing Microscope), Stuttgart 1953, p. 279.

## II. Basic Procedure for Examinations\*

Introductory Remarks: Check the centration of the source of light at frequent intervals with the Bertrand lens.

Check the position of axis  $A_2$  frequently with the auxiliary circles.

If the adjustments are disturbed by pronounced dispersion, use a good orange or red filter.

If necessary, use an eyepiece with iris diaphragm and cross-lines.

\* For the first exercises it is advisable to use rock sections orientated vertically to an optical symmetry axis, and to start accordingly with II, 5.

### 1. Finding the plane of symmetry.

- (a) Rotate on  $A_1$  until extinction.
- (b) Test whether on rotating on  $A_4$ , extinction is retained.  
Likewise for the other position of extinction.  
For the further operations, use the position with optimum extinction.
- (c) Alternating inclination on  $A_2$  and restoration of extinction through turning on  $A_1$  until extinction is retained on turning on  $A_4$ .
- (d) Take readings of the co-ordinates  $\alpha'_1$  and  $\alpha'_2$  ( $\alpha'_2$  is positive if read on the right auxiliary circle, and negative if read on the left).

### 2. Finding a second plane of symmetry.

- (a) Axes  $A_2$  and  $A_4$  in the normal position.
- (b) Index of  $A_1$  at  $\alpha'_1 - 90$ , if  $\alpha'_2$  positive  
"  $\alpha'_1 + 90$ , if  $\alpha'_2$  negative.
- (c) either: rotate  $A_1 \leftarrow$  to the left and lower  $A_2$  to the left until extinction is reached,  
or: rotate  $A_1 \rightarrow$  to the right and lower  $A_2$  to the right until extinction is reached. The smaller /  $\alpha'_2$  /, the smaller should be the rotations on  $A_1$ , and incline primarily on  $A_2$  only.  
Continued checking of the extinction by inclining on  $A_4$ .  
Take readings of the final co-ordinates  $\alpha''_1$  and  $\alpha''_2$ .

### 3. Determining the co-ordinates of a symmetry axis.

Either: with the aid of Wulff's net ruling.

- (a) Mark the North pole on an auxiliary sheet (tracing paper).
- (b) Mark the angles  $\alpha'_1$  and  $\alpha''_1$  on the outer circle in an anticlockwise direction from the North pole: points  $P'_1$  and  $P''_1$ .
- (c) Rotate the auxiliary sheet until  $P'_1$  coincides with the North pole of the net ruling; mark the angle  $\alpha'_2$  at the equator from the centre (to the right if  $\alpha'_2$  is positive, or to the left if  $\alpha'_2$  is negative): point  $P'_2$ .  
Repetition of the same procedure for the point  $P''_1$ ,  $\alpha''_2$  gives point  $P''_2$ .
- (d) Mark the point of intersection of the two great circles determined by  $P'_1, P'_2$  and  $P''_1, P''_2$ .
- (e) Rotate the auxiliary sheet to the right  $\rightarrow$  until S falls on the equator; read the angle of rotation from the North pole on the outer circle:  $\hat{\alpha}_1$ . Read the inclination of S towards the centre:  $\hat{\alpha}_2$  (positive on the right half, negative on the left).  
Or: by calculation with a slide rule.

$$(a) \tan (\alpha'_1 - \hat{\alpha}_1) = \frac{1}{\sin (\alpha'_1 - \alpha''_1)} \frac{\tan \alpha''_2}{\tan \alpha'_2} - \cotan (\alpha'_1 - \alpha''_1).$$

- (b) Thus giving  $\hat{\alpha}_1$ .

(c)  $\tan \alpha_2 = \frac{\tan \alpha'_2}{\cos (\alpha'_1 = \alpha_1)}$       The signs demand special attention!

$\alpha_1$  and  $\alpha_2$  are the co-ordinates of a symmetry axis.

**4. Setting a symmetry axis in the direction of observation.**

Set index  $A_1$  at  $\alpha_1$ , index  $A_2$  at  $\alpha_2$ , index  $A_3$  at 0.

Test: After rotating on  $A_3$  in each of the two positions of extinction, the extinction must be retained on inclining on  $A_1$ .

**5. Determining the characteristic angle of extinction.**

(a) Set axis  $A_1$  at 0.

(b) Rotate on  $A_3$  in the diagonal position.

(c) Determination of the character of the direction of oscillation SW-NE with gypsum red of the first order, with marked refraction with a quartz wedge or a compensator; note down the results.

(d) Rotate on  $A_3$  clockwise until first extinction; repeat four times; as accurately as possible: average of the readings  $\alpha_3$ .

(e) Set  $A_3$  at index  $\alpha_3 + 45^\circ$ .

(f) Set  $A_4$  at the value  $\pm 54.7^\circ$ .

With considerable differences in the light refraction of the mineral ( $n$ ) and that of the segment ( $n'$ ),  $A_4$  should be set to the value  $\frac{n}{n'} 54.7$ .

(g) Rotate on  $A_5$  anticlockwise until first extinction; repeat four times; as accurately as possible: average of the readings  $\alpha_5$ .

(h) Characteristic angle of extinction:

$\varphi =$  zero position of  $A_5 - \alpha_5$ .

**6. Determining the optical characteristic.**

(a) True angle of the optical axes: Determine  $\varphi$  in the standard diagram, read 2 V.

(b) Character of double refraction: upper sign in the diagram, if in 5c the character of the direction of oscillation was determined as being SW-NE: c (k), but the lower sign if the character determined was a (g).

(c) Orientate the 1st and 2nd centre lines in relation to the symmetry axis as set according to 4. above, as shown in the diagram. Accurate specification of the position of the centre lines in relation to the rock section normal by means of the angle  $\alpha_1$ , and in relation to edges, cleavages, and the border lines of twinings in the specimen by means of the angle  $\alpha_1$ .

**7. Repetition of the determinations (5 and 6 above) in a special diagram.**

(a) Selection of the test conditions  $\tilde{\varphi}$  and  $\gamma$  as functions of the angle  $\varphi$  in 5 h:

$\varphi$	$\bar{\varphi}$	$\chi$	$\varphi$	$\bar{\varphi}$	$\chi$
0 — 2	100	10	88 — 90	10	10
1 — 6	100	20	84 — 89	10	20
5 — 15	100	35	75 — 85	10	35
10 — 23	100	50	67 — 80	10	50
20 — 28	100	65	62 — 70	10	65
27 — 33	105	65	57 — 63	15	65
32 — 37	115	65	53 — 58	25	65
35 — 45	130	65	45 — 55	40	65

- (b) Set  $A_3$  at  $a_3 + \bar{\varphi}$  ( $a_3$ : average from 5 d).
- (c) Set  $A_4$  to the value  $\pm \chi$  or  $\pm \frac{n}{n'} \cdot \chi$ .
- (d) Rotate on  $A_5$  anticlockwise until first extinction; repeat four times; average  $a_3$ .
- (e) Extinction angle  $\varphi =$  zero position of  $A_3 - a_3$ .
- (f) Read 2 V from the special diagram. The cases in the same horizontal row of the table (a) above require the same diagram. The construction data for these are given on page 39 of the original paper.

Besides polished thin sections, individual mineral grains or crystals can be examined with the U-stage (see literature under Dimler-Stahmann).

With the method of H. Waldmann, the segments are replaced with a hollow glass sphere, in which grains of 1–11 mm. diameter can be measured morphologically and crystal-optically.

For examinations in incident light, the UTR 2 is available.

The following literature will facilitate introduction to the universal rotating stage methods, or give instructions for carrying out special examinations.

- Berek, M. Mikroskopische Mineralbestimmung mit Hilfe der Universaldrehtischmethoden (Berlin 1924) Gebr. Borntraeger.
- Reinhard, M. Universaldrehtischmethoden (Basel 1931), Wepf & Cie, (mit besonderer Berücksichtigung der Plagioklasbestimmung).
- Nikitin, W. Die Fedorow-Methode. (Berlin 1936).
- Emmons, R. C. The Universal Stage (with five axes of rotation). Mem. Geol. Soc. Amer. 8, 1943.  
Ausführliche Darstellung der I-T-Variationsmethode.
- Chudoba, K. Die Feldspäte und ihre praktische Bestimmung. (Stuttgart 1932).

- Nikitin, W. La méthode universelle de Fédorow.  
Traduction française par L. Duparc et V. de Dervies.  
(Genève, Paris, Liège 1914).
- 
- Chomard, L. Theorie et pratique de la méthode de Fédorow.  
Procédé classique et méthode analytique générale.  
Ann. Mines (13), Memofres 5, (1934), P. 153-218.
- Berek, M. Grundsätzliches zur Bestimmung der optischen Indikatrix mit Hilfe  
des Universaldrehtisches.  
Schweiz. Mineral.-Petrogr. Mitt. 29, (1949), S. 1-18.
- Emmons, R. C. A modified universal stage. Amer. Mineral. 14, (1943), pp. 441-461.
- Dimler, R. J. und  
Stahmann, M. A. A mount for the universal stage study of fragile materials,  
Amer. Mineral. 25, (1940), pp. 502-504.
- Waldmann, H. Glashohlkugel für Kristall- und Edelsteinuntersuchungen.  
Schweiz. Mineral.-Petrogr. Mitt. 27, (1947), S. 472-520.
- Tröger, E. Nomogramm zur Reduktion von Kippwinkeln am Universaldrehtisch.  
Cbl. Mineral. A, (1939), S. 472-520.
- Chudoba, K. Bestimmung der Plagioklase in Doppelzwillingen nach dem Albit- und  
Karlsbadergesetz bei Untersuchungen mit Hilfe des Universaldreh-  
tisches. N. Jb. Mineral. A. Beil. - Bd. 63, (1932), S. 267-278.
- Tertsch, H. Die Ermittlung der kristallogr. Orientierung der Plagioklasoptik ohne  
Verwendung von Spaltrissen oder Flächenspuren.  
Z. Kristallogr. 105, (1944), S. 330-334.
- Köhler, Al. Drehtischmessungen an Plagioklaszwillingen von Tief- und Hoch-  
temperaturoptik. Mineral.-Petrogr. Mitt. 53, (1942), S. 159-179.
- Rittmann, A. Die Zonenmethode. Ein Beitrag zur Methodik der Plagioklasbestim-  
mung mit Hilfe des Théodolitischen.  
Schweiz. Mineral.-Petrogr. Mitt. 9, (1929), S. 1-46.
- Bäier, E. Der Zwillingsfall als Spezialfall zentrischer Verdrehung zweier Kristalle.  
Z. Kristallogr. 87, (1934), S. 306-325.
- Nikitin, W. Korrekturen und Vervollständigungen der Diagramme zur Bestim-  
mung der Feldspäte nach Fedorows Methode.  
Mineral.-Petrogr. Mitt. (Leipzig) 44, (1933), S. 117-160.
- Melon, J. Essai de détermination des propriétés optiques d'un minéral par la  
mesure, en lumière parallèle oblique, des retards en différents points  
d'une lame cristalline. Ann. Soc. Geol. Belgique (Liège) 57, M 3-M 108,  
(1933-1934).

- Mosebach, R. Über die Bestimmung der Brechzahlen doppelbrechender Minerale im gedeckten Dünnschliff, Nachr. Akad. Wiss. Göttingen, Math.-Phys. Kl. (1947), S. 20–22.
- Burri, Conrad Das Polarisationsmikroskop. – Verlag Birkhäuser, Basel (1950).
- Nieuvenkamp, W. M. Berek's Methode der charakteristischen Gangunterschiedsverhältnisse, Proc. k. Akad. Wetensch. Amsterdam **30**, (1927), S. 534–542.
- Dodge, T. A. The determination of optic angle with the universal stage. Amer. Mineral. **19**, (1934), pp. 62–75.
- Schumann, H. Erweiterung der konoskopischen Beobachtungsweise durch den Drehtisch. Fortschr. Mineral. **21**, (1937), S. 102–105.
- Schumann, H. Über den Anwendungsbereich der konoskopischen Methodik. *ibid* **25**, (1941), S. 217–252.
- Schumann, H. Demonstration allseitig drehbarer konoskopischer Interferenzbilder mittels Mikroprojektionsapparates und Drehtisches. Cbl. Mineral. Abt. A (1938), S. 120–127.
- Sander, Bruno Gefügekunde der Gesteine (1930). Wien.
- Dietzel, A. Die Identifizierung der Entglasungsprodukte von techn. Natron, Kalk, Kieselsäure, Gläsern. Glastechn. Ber. **9**, (1931), S. 307–309.
- Doeglas, Th. J. A rapid microscopic method for distinguishing quartz from untwinned oligoclase-andesine. Amer. Mineral. **21**, (1936), pp. 531–532; **25**, (1940), pp. 288–296.
- Duparc, L. und Reinhard, M. La détermination des plagioclases dans les coupes minces. Genève 1924.
- Plato, W. Bestimmung kristallogr. Daten kristallis. geschliffener Edelsteine. Z. wiss. Mikrosk. **57**, (1940/41), S. 137.
- Drescher, F. K. Ein neues großes Universal-Drehinstrument. Fortschr. Mineral. **15**, (1931), S. 319–320.
- Ehrenberg, H. Ein neuer Universaldrehtisch für Anschliffe und seine Anwendungsmöglichkeiten. Erzmetall **3**, (1950), S. 65–69.

**ERNST LEITZ GMBH WETZLAR**  
Germany

Branch Works: Ernst Leitz (Canada) Ltd., Midland, Ontario.

55-9/Engl.

IV/56/AX

Printed in Germany  
Scharfdruck Wetzlar