

VICKERS
FIFTY-FIVE
MICROSCOPE
by COOKE



M 550002

With the addition of
M 550013 high power
tungsten filament lamp


VICKERS FIFTY-FIVE MICROSCOPE

INDEX

	<i>Page</i>
Summary of instrument	4
Sectional drawings	8
General description	12
Automatic integrating photographic timer	14
35 mm. Camera unit	14
Automatic magnification indicator	16
Eyepiece magnification changers	16
Multiple objective carriers	16
Incident and transmitted illumination	18
Phase contrast equipment	20
Polarizing equipment	22
Dark ground equipment	24
Oblique illumination unit	25
Macro low power equipment	26
Micro hardness testing equipment	28
Photographic equipment	29
Electrical equipment	30
List of objectives	31
List of eyepieces	31
Photographic accessories	33
General accessories	34

COOKE TROUGHTON & SIMMS LTD.
HAXBY ROAD
YORK • ENGLAND





VICKERS FIFTY-FIVE MICROSCOPE

M 550001



M550001 Vickers Fifty-Five Microscope with gliding stage (joy stick control), normal incident illuminator, monocular eyepiece, sextuple carrier, magnification changer, xenon lamp (interchangeable with high pressure mercury vapour lamp), electrical supply cabinet (wall mounting), focusing screen and double plate holder for $\frac{1}{2}$ plate with two $\frac{1}{4}$ plate adaptors and hood, adjustable lamp pillar carrying xenon or mercury vapour lamp, centring controls to lamp, lamp condenser on focusing mount, water trough, colour filter, neutral filter and diffusing screen. Dimensions: Height $30\frac{1}{2}$ in., Width 44 in., Depth 24 in., Weight: 177 lbs.

COOKE TROUGHTON & SIMMS LTD.




VICKERS FIFTY-FIVE MICROSCOPE

M 550002




M550002 Vickers Fifty-Five Microscope with centring, traversing, rotating and gliding stage, normal incident illuminator, binocular eyepiece, sextuple carrier, magnification changer with analyser unit, xenon lamp (interchangeable with high pressure mercury vapour lamp), transmitted light bracket, motorised 35 mm. camera unit, integrating photographic timer, motorised focal plane shutter, wood cabinet with built-in electrical supply, focusing screen and double plate holder for $\frac{1}{2}$ plate with two $\frac{1}{4}$ plate adaptors and hood, adjustable lamp pillar carrying xenon or mercury vapour lamp, centring controls to lamp, lamp condenser on focusing mount, water trough, colour filter, neutral filter and diffusing screen. Dimensions: Overall Height 61 in., Cabinet Floor Space 44 in. \times 24 in. Weight: Microscope 204 lbs., Cabinet 186 lbs.

COOKE TROUGHTON & SIMMS LTD.



- | | |
|-------------------------------|---|
| 1. Stage Clips | 12. Rack and Pinion (Coarse Motion) |
| 2. Specimen | 13. Zoom Projection Eyepiece |
| 3. Gliding Stage | 14. Mirrors (2) |
| 4. Micrometer Stage | 15. Stage Traverse Micrometers |
| 5. Stage Support | 16. Binocular Eyepiece |
| 6. Multiple Objective Carrier | 17. Selector Switch (Visual only—35 mm. photo and visual—Macro) |
| 7. Illumination Box | 18. Focal Plane Shutter |
| 8. Slow Motion Carriage | 19. Selector Rod (Zoom Eyepiece—Macro—35mm. Corrector Lens) |
| 9. Magnification Changer | 20. Focusing Screen |
| 10. Weight Relieving Spring | 21. 35 mm. Camera |
| 11. Slow Motion Transfer Gear | 22. Anti-Vibration Mountings |



- | | |
|---|--|
| 23. Stage Clamp Screw | 34. Magnification Changer Scale |
| 24. Field Iris Control | 35. Weight Relieving Spring for Xenon Lamp |
| 25. Photomultiplier Cell | 36. Xenon Lamp |
| 26. Microscope Block | 37. Reflector Unit |
| 27. Subsidiary Coarse Motion Pinion Control | 38. Condenser Focusing Control |
| 28. Fine Motion Heads | 39. Condenser Bracket |
| 29. Clamp for Coarse Motion | 40. Lamp Centring Control Heads |
| 30. Coarse Motion Milled Head | 41. Clamp Lever for Lamp Slideway |
| 31. Pick-off Prism for Photomultiplier | 42. Xenon Lamp Casing |
| 32. Objective Setting Scale | 43. Lamp Mains Switch |
| 33. Final Magnification Scale | 44. Multi-pin Plug |

**NEW
and
SUPERIOR FEATURES**

Automatic Integrating Photographic Timer

Zoom Projection Eyepiece

Eyepiece Magnification Changer

35 mm. Camera, Manual or Motorised Control

Multiple Objective Carriers

Automatic Magnification Indicator

Built-in Polarizing Unit

Simultaneous Incident and Transmitted Illumination

Pneumatic Micro-Hardness Testing Equipment

Incident Phase Contrast Unit

Instrument mounted on Anti-Vibration Pads


All Illumination Techniques applicable without disturbing specimen



M550300 Micrometer and Gliding Stage



M551515 Gliding Stage



M550017
Automatic Integrating
Photographic Timer

Wall Mounting
Power Supply Unit



OPE

This break-down display of major components can serve as a guide when ordering an instrument to suit individual requirements.

M550012 Eyepiece Magnification Changer
with built-in Analyser Unit




M550011
Eyepiece Magnification Changer




M550680 Binocular Body



M550670
Monocular Body



M550015 35 mm. Camera Unit



M551850 Focal Plane Shutter



VICKERS FIFTY-FIVE MICROSCOPE

GENERAL DESCRIPTION

The Vickers Fifty-Five Microscope is of the inverted type, the specimen being positioned above the objectives. It is a bench type instrument, with built-in anti-vibration mountings, and versions can be supplied either to fit an existing bench installation or complete with cabinet containing electrical equipment.

Exceptional stability and freedom from vibration under the highest magnifications is achieved due to the manner in which the objectives and mount, together with the fine motion mechanism, are connected to the stage support.

The standard object stage has a gliding top plate which responds smoothly to the lever control. For precise measurement and polarizing work a more comprehensive stage is available, having, in addition to the gliding motion, micrometer-controlled lateral and transverse movements, together with a graduated rotary movement so arranged as to be always centred to the optical axis.

An important feature is that once the object under examination has been placed in position on the stage it need not be moved when changing the type of illumination technique employed to view it; micro hardness testing can also be performed in situ.

The relationship between the angle of the eyepiece and the viewing screen has been so designed that the observer may rapidly alternate between this and the eyepiece when desired.

The main source of illumination is a High Pressure Xenon Lamp which is interchangeable on a vee slide with a Mercury Vapour Lamp. The lamp and lamp housing are designed to move vertically under the control of a spring-loaded device so that the light may be fed into any of the various illuminators available. Full control over the illumination beam is ensured by the provision of centring adjustments, iris diaphragm, focusing lamp condenser and built-in filters. When required, mixed illumination is available, the xenon source providing normal incident light and a tungsten filament lamp supplying transmitted light. The tungsten filament lamp is suitable for most transmitted light work.

For incident illumination, light enters the side of an incident illuminator unit via a field iris, and is reflected upwards through the selected objective, mounted on a rotatable objective carrier, to the specimen. The image-forming rays then return via a magnification changer having alternative settings providing magnifications of $1.0\times$, $1.4\times$ and $2.0\times$. The light path may now be deflected into a monocular or binocular head carrying conventional eyepieces and/or be allowed to continue to the photographic focusing screen. The choice of light path is governed by a selector switch, on the front of the focusing block, which may be set to one of three positions, "Visual only", "35 mm. Photo and Visual", or "Macro-Photo".

In the "Visual only" position all the light is directed to the viewing eyepiece. The magnification may be calculated as a product of the objective, magnification changer and eyepiece powers. When the binocular head is in use, tube-length compensators are fitted to each eyepiece tube to allow correction for varying interocular separations.

When the selector switch is in the "Photo and Visual" position, 90 per cent of the available light is passed to the photographic screen, and the remainder deflected to the visual head. This is to enable accurate focusing and centring of the image. Photographs on half or quarter plates ($7'' \times 5''$ and $5'' \times 4''$ respectively in North America) may then be produced by substituting sensitive material in place of the viewing screen.

The third position of the selector switch "Macro-Photo" enables the light beam to be projected directly on to the screen allowing macro examination and photography, either by means of transmitted or incident illumination, macro lenses and illuminators being substituted for the magnification changer and micro-objective assembly.

The swing-out mirror bracket is necessary when using transmitted light techniques.

The instrument is particularly well equipped for work in polarized light, since all the necessary polarizing components can be built into a special magnification changer.

The Vickers Fifty-Five Microscope incorporates a panoramic Projection Eyepiece which, in conjunction with the magnification changer and a series of objectives, allows a comprehensive range of screen magnifications to be produced, and the magnification value may be

read directly on the Zoom control knob, which is calibrated with three scales, the outer scale is set to indicate the power of the objective in use, the inner scale to the value selected on the magnification changer, while the centre scale, a graduated arc, gives a direct reading of the final screen magnification as the control is rotated against it.

The range covered by the "Zoom" projection eyepiece is from 24 diameters with a 3.5 \times objective, and the magnification changer set at 1, to 2,800 diameters with a 140 \times objective and the magnification changer set at 2.0. A.S.T.M. recommended magnifications are shown in red.

The objectives used for incident light, transmitted light and dark field microscopy are mounted on sextuple or quintuple carriers, which are quickly interchangeable. All incident light objectives used on this instrument are computed to work at infinite tube length. 160 mm. objectives are used for transmitted light work, and their quintuple mount incorporates a corrector lens. All objectives except 3.5 \times and 6.0 \times are par-central and par-focal.

When dark field optics are used the necessary patch-stop is incorporated on a slide mounted on the front of the lamp condenser. This slide also carries an annulus which may be matched to phase rings for incident phase contrast, matching being carried out by the use of a focusing Bertrand lens fitted into the magnification changer. The incident phase contrast unit complete with the illumination box replaces the normal incident illuminator and provision is then made for positive and negative phase contrast, dark ground illumination and normal incident illumination. For transmitted phase contrast a condenser annulus changer together with 160 mm. tube length phase objective is required.

A Photographic Shutter can be incorporated into the instrument and may be used for a range of exposure times when using either a plate camera or a 35 mm. camera. When the 35 mm. camera is inserted into the instrument it automatically positions a selector switch which places a corrector lens into the light path. On removal of the 35 mm. camera the selector may be set to either the projection eyepiece or "Macro" position.


The Automatic Integrating Photographic Timer removes the need for trial exposures. A photo-cell, fed from the optical path, integrates the image-forming light and measures its intensity, terminating the exposure when a pre-determined amount of light energy has passed through. When the switches for the film speed and photometer scale have been set correctly, depression of the spring switch marked "Exposure" controls the photographic shutter and automatically gives the correct exposure to the sensitive photographic material. The Exposure Duration Indicator gives a visual warning that the shutter is open and an exposure is in progress.

AUTOMATIC INTEGRATING PHOTOGRAPHIC TIMER

The Automatic Integrating Timer (developed in collaboration with Dr. R. Barer and Mr. J. Underwood, Department of Human Anatomy, Oxford University) will accommodate film speeds over the range 5 A.S.A. to 3,200 A.S.A. and automatically timed exposures can be obtained at shutter speeds from 1/10th second to several minutes. Shutter speeds up to 1/500th second can be set on the shutter, the exposure time having been calculated from the reading of the built-in photometer. The shutter is automatically rewound after each exposure, and if the 35 mm. camera is in use, the film in the camera is transported one frame.


It is possible to follow the progress of an exposure visually by means of the exposure duration indicator.

The Automatic Timer can also be used with the macro equipment.





M551850 Focal plane shutter



M550015 35 mm. Camera unit showing cassette



- M550017** Automatic Integrating Photographic Timer (including focal plane shutter).
- M551850** Focal Plane Shutter.
- M550015** 35 mm. Camera Unit (Motorised).
- M550014** 35 mm. Camera Unit (Manual Control only).
- M551854** Additional Laboratory Cassette (with film advance mechanism).

When ordering it is essential to give particulars of electricity supply.

AUTOMATIC MAGNIFICATION INDICATOR



Automatic magnification indicator allows direct reading from 24 to 2,800 diameters. A.S.T.M. recommended magnifications are shown in red.

EYEPIECE MAGNIFICATION CHANGERS




The magnification changer contains a rotating lens system giving 1.0 \times , 1.4 \times and 2.0 \times magnifications, consequently a large choice of eyepieces is unnecessary.

- M550011 Eyepiece Magnification Changer.
- M550012 Eyepiece Magnification Changer with Built-in Analyser Unit.

MULTIPLE OBJECTIVE CARRIERS




- M550560 Sextuple carrier and mount for infinite tube length normal incident objectives (supplied with the instrument).
- M550550 Quintuple carrier and mount for infinite tube length 4 mm., 8 mm. and 16 mm. dark ground objectives with Catoptric Condensers and two normal incident objectives.
- M551545 Quintuple carrier and mount with built-in 2 \times corrector lens for 160 mm. tube length transmitted light objectives.



Automatic magnification indicator



M550011



M550012




M550560



M550550




M551545



INCIDENT AND TRANSMITTED ILLUMINATION

Incident light arrangement for the examination of a metallurgical specimen.




EQUIPMENT FOR TRANSMITTED LIGHT

- M550500 Substage Swing-out Bracket.
- M001376 Centring Condenser Mount.
- M001382 Abbe Condenser (2 lens).
- M001383 Aplanatic Condenser (4 lens).
- M001386 Achromatic Condenser N.A. 1.00
- M001391 Achromatic Oil Immersion Condenser N.A. 1.30.


Transmitted light arrangement, with substage swing-out bracket, for the examination of a biological specimen.

SIMULTANEOUS INCIDENT AND TRANSMITTED ILLUMINATION




Simultaneous incident and transmitted illumination can be achieved by the addition of a high power tungsten filament lamp.

M550013 High Power Tungsten Filament Lamp.



PHASE CONTRAST EQUIPMENT FOR INCIDENT LIGHT



Incident phase contrast complete with illumination box allowing positive and negative phase contrast, dark ground and normal incident illumination.

M550570 Incident Phase Unit.


OBJECTIVES FOR INCIDENT LIGHT PHASE CONTRAST

These objectives are corrected for use on uncovered specimens and infinite tube length.

Achromatic Objectives	Power	Numerical Aperture
M022352	15×	0.25
M022452	30×	0.5
M022952	50×	0.8
M022652	140×	1.3

A complete list of eyepieces appears on page 31

PHASE CONTRAST EQUIPMENT FOR TRANSMITTED LIGHT



Transmitted phase contrast showing substage swing-out bracket and phase condenser annulus.

- M550500 Substage swing-out bracket.
- M555625 Phase contrast unit complete with condenser and four annular diaphragms.
- M555642 Phase contrast unit as M555625 but with long working distance condenser.
- M551545 Quintuple objective carrier with 2 × corrector lens.

OBJECTIVES FOR TRANSMITTED LIGHT PHASE CONTRAST

These objectives are corrected for use with cover glass and 160 mm. tube length.

Achromatic Objectives	Power*	Numerical Aperture
M022205	20 ×	0.25
M022405	40 ×	0.5
M022505	80 ×	0.65
M022605	190 ×	1.3
Fluorite Objective		
M023605	90 ×	0.95

*Including corrector lens.

A complete list of eyepieces appears on page 31

POLARIZING EQUIPMENT FOR INCIDENT AND TRANSMITTED LIGHT

The monocular eyepiece with swing-out focusing Bertrand lens and a magnification changer with built-in analyser are required for work with polarized light. The polaroid analyser and quartz sensitive tint plate incorporated in the magnification changer rotate together, and measurement of this rotation can be obtained from a graduated scale with a range 0 to 105 degrees which can be read, against a vernier, to 12 minutes. The polaroid analyser and the quartz sensitive tint plate can be withdrawn from the light path independently.

A unique feature of the magnification changer unit is the addition of a rotating slot designed to take a range of compensators, the rotary movement of 360° being divided in degrees, and read by a vernier to 6 minutes.


INCIDENT LIGHT

For incident light work a graduated polarizing cap, which fits on the incident illuminator, is required.

TRANSMITTED LIGHT


Additional equipment necessary for transmitted light consists of the quintuple transmitted light objective carrier, fitted with a $2\times$ corrector lens for 160 mm. tube length objectives, and a substage swing-out mirror bracket with a condenser and polarizing substage.

- M551897 Monocular eyepiece with Bertrand lens.
- M550012 Magnification changer with analyser unit.
- M551545 Quintuple objective carrier with $2\times$ corrector lens.
- M551040 Polarizing substage attachment.
- M550500 Substage swing-out bracket.
- M552065 Graduated polarizing cap (for incident light only).



Instrument set up for incident polarized work

Arrangement for transmitted polarized work



ACCESSORIES for incident and transmitted work.

- M552075 Quartz wedge, 6 orders (non-graduated)
- M552085 Mica $\frac{1}{4}$ wave plate.
- M552090 de Sénarmont compensator.
- M552080 Quartz wedge, 6 orders (graduated).
- M552056 Nakamura half-shadow plate.
- M550950 Compensator plate (elliptic).

CONDENSERS for transmitted polarizing work.

- M007882 Abbe condenser (2 lens).
- M007884 Aplanatic condenser (4 lens).
- M007886 Achromatic condenser N.A. 1.0.
- M007891 Achromatic oil immersion condenser N.A. 1.3.


ACHROMATIC OBJECTIVES FOR POLARIZING WORK

- A** Objectives corrected for use with uncovered specimens and for infinite tube length.
- B** Objectives corrected for use with specimens having a cover glass 0.007 in. (0.18 mm.) thick and for 160 mm. tube length. The quintuple objective carrier incorporates a 2 \times corrector lens.

A	B	Type	Numerical Aperture	Primary Magnification		Working Distance mm.	
				A	B*	A	B
M023054	—	Achro.	0.05	3.5 \times	—	69	—
M022054	M022003	„	0.10	6 \times	6 \times	38	43
M022154	M022103	„	0.15	10 \times	10 \times	14	17
M022354	M022303	„	0.25	15 \times	20 \times	5	5
M022454	M022403	„	0.50	30 \times	40 \times	1.42	1.52
—	M022503	„	0.65	—	80 \times	—	0.71
—	M022903	„	0.85	—	80 \times	—	0.43
M022954	—	„	0.80	50 \times	—	0.48	—
M023154	—	„	0.85	85 \times	—	0.28	—
M022654	M022603	Achro. oil	1.30	140 \times	190 \times	0.41	0.17

*Including corrector lens.

A complete list of eyepieces appears on page 31



EQUIPMENT FOR INCIDENT DARK GROUND ILLUMINATION

For dark ground work a quintuple carrier and mount to take the infinity corrected 4 mm., 8 mm., and 16 mm. dark field objectives with catoptric condenser, is required. This carrier also accommodates two normal incident objectives. The patch stop is situated on a slide in front of the lamp condenser.

- M550550** Quintuple carrier for dark ground objectives.
- M551896** Catoptric condenser for use with M023452 and M024152 dark ground objectives.
- M551895** Catoptric condenser for use with M023252 dark ground objective.

Equipment for transmitted dark ground illumination is listed on page 33

Achromatic Objectives for Dark Ground Illumination

Achromatic Objectives	Type	Numerical Aperture	Primary Magnification	Working Distance mms.
M023252	Dark Ground	0.25	15×	5
M023452	„ „	0.50	30×	1.42
M024152	„ „	0.65	50×	0.75


M023252 dark ground objective without the patch stop will serve the same purpose as the M022352 15× achromatic objective.

M023452 and M024152 may also be used for bright field work.

A complete list of eyepieces appears on page 31

OBLIQUE ILLUMINATION UNIT

Incident light arrangement with oblique illumination unit.





The oblique illumination unit slides into the filter holder situated in front of the lamp condenser. To regulate the obliquity of the illumination beam, provision is made for the aperture diaphragm to be decentred and rotated. Scales are provided on both the decentring and rotating movements to enable correct and repetitive settings of the oblique beam.

This type of illumination is designed to give exceptional contrast and to increase depth perception.


M551960. Oblique Illumination Unit.

Illustration of the viewing hood which is supplied with the instrument.





*Incident light arrangement
for 5× macro objective.*



*Incident light arrangement
for 10× and 15× macro
objectives.*

MACRO

LOW POWER EQUIPMENT



For both incident and transmitted macro examination the eyepiece magnification changer is removed and replaced by the macro base unit which carries the 5×, 10× and 15× objectives, and a ground glass screen is inserted into the lamp condenser mounting.

MACRO incident illumination

For incident illumination two incident illuminators are offered, one for use with the 10× and 15× objectives, the other, together with a supplementary projector lens, for use with the 5× objective.


An oblique illumination mirror is also available.

MACRO transmitted illumination


For transmitted illumination the substage swing-out mirror bracket and auxiliary condenser, together with a projection lens, are required.

MACRO OBJECTIVES

	Focal Length (inches)	Aperture	Magnification
M551325	5.0	f/4.5	5×
M551340	2.6	f/4.5	10×
M551360	1.9	f/3.5	15×



Oblique incident illumination mirror.



Arrangement for transmitted macro work.

COMBINATIONS OF MACRO EQUIPMENT

Macro Objectives and types of Illumination	M551270 Macro Base Unit	M551795 10× & 15× Incident Illuminator	M551815 5× Incident Illuminator	M551310 Oblique Illuminator Mirror	Ground Glass Screen	M551825 Auxiliary Condenser, Transmitted	M551990 Projection Lens for Transmitted	M551830 Projection Lens for 5× Incident	M550500 Swing-out Mirror Bracket
15× Incident	✓	✓			✓				
15× Transmitted	✓				✓	✓	✓		✓
10× Incident	✓	✓			✓				
10× Transmitted	✓				✓	✓	✓		✓
5× Incident	✓		✓		✓			✓	
5× Transmitted	✓				✓	✓	✓		✓
15× Inc. Oblique	✓			✓					
10× Inc. Oblique	✓			✓					
5× Inc. Oblique	✓			✓					
10× & 15× Inc.	✓	✓			✓				
10× & 15× Trans.	✓				✓	✓	✓		✓
5×, 10× & 15× Inc.	✓	✓	✓		✓			✓	
5×, 10× & 15× Trans.	✓				✓	✓	✓		✓


MICRO HARDNESS TESTING EQUIPMENT

Micro hardness testing can be efficiently carried out by non-specialised personnel, the apparatus having been designed to operate automatically under predetermined loads between 5 gms. and 200 gms. Specimens may be tested without the use of special holders or chucks.

The micro hardness testing equipment consists of a combined indenter and microscope objective in which the housing for the optical lens system is fixed, and only the diamond indenter, mounted on a thin rubber diaphragm, moves.


After the selection of the required load, which is read against the index on the cylinder of the transmitter, a lever is depressed, and the correct load at a constant rate for that load is automatically applied to the indenter by pneumatic pressure. When the lever is allowed to return, the pneumatic pressure is released, and the indenter is restored to its original position. The indenting operations are entirely impersonal, and variations in the rate of operation of the lever, which may be expected between different operators, do not influence the final results.

M550016 Micro Hardness Testing Equipment, including micro hardness transmitter, microscope indenter objective, 4 mm. achromatic objective and centring eyepiece, with filar micrometer.




Arrangement for micro hardness testing.

Micro hardness transmitter, microscope indenter objective and filar micrometer eyepiece.





M552095



M552100

PHOTOGRAPHIC EQUIPMENT

- M552095 M.P.P. camera back (to take 5 in. \times 4 in. plates), wood adaptor, and focusing hood.
- M552100 M.P.P. camera back, wood adaptor, focusing hood, and polaroid land (5 in. \times 4 in.) film holder.
- M552105 M.P.P. camera back, wood adaptor, focusing hood, and polaroid roll film back.

Photographic accessories are listed on page 33.

M552105



ELECTRICAL EQUIPMENT

The power supply unit is available in two forms; one designed for wall mounting, the other built into the wood cabinet.


The xenon lamp, which is the normal source of illumination, is fitted with a control and starter unit suitable for 50 or 60 cycles and 110, 200 and 240 mains voltage. A high pressure mercury vapour lamp with choke can be supplied as an alternative to the xenon lamp, and is interchangeable on a V-slide attachment. The mercury vapour lamp is only suitable for 220 to 240 mains voltage. Either the xenon or the mercury vapour lamp can be used for transmitted or incident light techniques.

With the addition of a high power tungsten filament lamp (48 watt), operating in conjunction with either the xenon or mercury vapour lamp, simultaneous incident and transmitted light can be obtained. The filament lamp is supplied with a control panel, mains switch, rheostat and transformer, which can be built into the instrument.

- M550013** High Power Tungsten Filament Lamp.
- M550018** Mercury Vapour Lamp (includes choke).
- M550790** Xenon Lamp (included with microscope).

When ordering it is essential to give particulars of electricity supply.

When electrical equipment not of our manufacture is supplied, our liability, in respect of any defect in or failure of the articles supplied, or for any consequential loss, injury, or damage, is limited to the benefit of any guarantee, condition or warranty given to us by the supplier or manufacturer, and then only to the extent to which we can enforce the same.



20-way plug connecting the instrument to the built-in power supply unit.

COOKE ACHROMATIC, APOCHROMATIC AND FLUORITE OBJECTIVES

A—Objectives for uncovered specimen. Infinite tube length.

B—Objectives for covered specimens, 160 mm. tube length used with transmitted light nosepiece containing 2× corrector lens.

		Strain free objectives for polarized light		Type	Normal Focal Length mm.	N.A.	Primary magnification		Working Distance mm.	
A	B	A	B				A	B*	A	B
M023052		M023054		Achro.	68	0.05	3.5×		69	
M022052	M022001	M022054	M022003	„	33	0.10	6×	6×	38	43
M022152	M022101	M022154	M022103	„	25	0.15	10×	10×	14	17
M022352	M022301	M022354	M022303	„	16	0.25	15×	20×	5	5
M022452	M022401	M022454	M022403	„	8	0.5	30×	40×	1.42	1.52
—	M022501	—	M022503	„	4	0.65	—	80×	—	0.71
—	M022901	—	M022903	„	4	0.85	—	80×	—	0.43
M022952	—	M022954	—	„	4	0.8	50×	—	0.48	—
M023152	—	M023154	—	„	3	0.85	85×	—	0.28	—
M022652	M022601	M022654	M022603	Achro. Oil	2	1.3	140×	190×	0.41	0.17
M024052	M024001	—	—	Apo.	16	0.3	17×	20×	39	5
M023852	M023801	—	—	„	4	0.95	50×	80×	0.3	†
M023352	—	—	—	„	3	0.95	85×	—	0.18	—
M023752	M023701	—	—	Apo. Oil	2.2	1.32	115×	160×	0.3	0.12
M023552	M023501	—	—	Fluorite Oil	1.8	1.3	140×	190×	0.3	0.12
M023252	—	—	—	Achro. Dark Ground	16	0.25	15×	—	5	—
M023452	—	—	—	„	8	0.5	30×	—	1.42	—
M024152	—	—	—	„	4	0.65	50×	—	0.76	—

Transmitted light phase contrast

—	M022205	—	—	Achro.	16	0.25	—	20×	—	5
—	M022405	—	—	„	8	0.5	—	40×	—	1.52
—	M022505	—	—	„	4	0.65	—	80×	—	0.71
—	M022605	—	—	Achro. Oil	2	1.3	—	190×	—	0.17
—	M023605	—	—	Fluorite Oil	3.75	0.95	—	90×	—	0.23

* Including corrector lens. † With correction collar for cover glasses.


Any 160 mm. tube length objective with a 34 mm. shoulder length can be used with the transmitted light nosepiece.

M550560 Sextuple carrier and mount for infinite tube length normal incident objectives is supplied with the instrument.


EYEPIECES

HUYGHENS			COMPENSATING			KELLNER		
Single	Paired	Power	Single	Paired	Power	Single	Paired	Power
M001506	M001507	6×	*M001526	*M001527	6×	M001541	M001542	10×
M001511	M001512	8×	*M001531	*M001532	8×	M001546	M001547	15×
M001516	M001517	10×	*M019536	*M019537	10×			
			*M001551	*M001552	20×			

* Can be used with spectacles.



*Diatom
craspedodiscus
coscinodiscus
Ehrenberg*



*Flow lines in section
of forged chain link*



Copper—copper phosphide eutectic

PHOTOGRAPHIC ACCESSORIES

- M505705 Photographic shutter.
- M505890 Screen focusing magnifier.
- M505664 Double plate holder for 7 in. \times 5 in. plates.
- M505689 Adaptors (2) for 5 in. \times 4 in. plates for M505664.
- M505647 Double plate holder for $\frac{1}{2}$ plates.
- M505960 Adaptors (2) for $\frac{1}{4}$ plates for M505647.
- M505790 Adaptors (2) for 12 cm. \times 9 cm. plates for M505647.
- M505931 Film holder for M505664.
- M505932 Film holder for M505689.
- M505933 Film holder for M505647.
- M505934 Film holder for M505960.
- M505935 Film holder for M505790.
- M551610 Frame for focusing screens.
- M505904 Focusing screen, ruled with cross lines.
- M505906 Focusing screen, with horizontal and vertical lines intersecting at centre of plate, each graduated 100 mm.
- M505907 Focusing screen, with grid 10 cm. \times 10 cm. ruled in mm., with break at intersection of 1 cm. lines.
- M505908 Focusing screen, with grid 10 cm. \times 10 cm. ruled in cm.

Other ruled focusing screens for measurement of particles, determination of grain size, etc. (Porton Globe and Circle, McQuaid and Ehn hexagon grain size screens), are available. Particulars on request.

It should be noted that although plate holders M505664 and M505931 take 7 in. \times 5 in. plates and film, the size of the picture will be less than this.

EQUIPMENT FOR TRANSMITTED DARK GROUND ILLUMINATION

- M001396 Dark ground condenser for transmitted light.
- M001397 Funnel stop for objectives of N.A. above 1.00 when used with M001396.
- M001394 Low power dark ground condenser for transmitted light fixed focus, 0.7—0.8 N.A.

ACCESSORIES

- M550670 Monocular body
- M550680 Binocular body
- M551515 Gliding stage.
- M550300 Micrometer and gliding stage.
- M001195 Small bottle of non-drying immersion oil, as supplied with each immersion objective, (ALP₁) R.I. 1·524.
- M001196 8 oz. bottle of non-drying immersion oil (ALP₁) R.I. 1·524.
- M001197 Small bottle of non-drying immersion oil (ALP₂) R.I. 1·515.
- M001198 8 oz. bottle of non-drying immersion oil (ALP₂) R.I. 1·515.
- M001581 Eyepiece micrometer of glass, 1 cm. divided into 100 parts.
- M001586 Stage micrometer of glass, 1 mm. divided into 100 parts.
- M001591 Stage micrometer of metal, 1 mm. divided into 100 parts.
- M001555 Filar micrometer 10× eyepiece with travelling web reading against a millimeter scale, and read directly to 0·01 mm. on a micrometer drum with estimation to 0·002 mm.
- M001556 Filar micrometer as M001555, but with 15× eyepiece.
- M551440 Cabinet and fittings.

SPARES

- M505970 Water trough.
- M505608 Heat absorbing filter.
- M505733 Xenon bulb X.B.O. 162 Neron Ltd.
- M551649 Mercury vapour bulb H.B.O. 200 Neron Ltd.
- E855 Tungsten filament bulb 6 volts, 8 amps (Siemens).
- M505795 Soft plastic dust cover (supplied with instrument).

Details of surface topography techniques will be provided on request.

As designs are constantly subject to revision the particulars listed throughout, may not be final in all respects.

COOKE TROUGHTON & SIMMS LTD.
A MEMBER OF THE VICKERS GROUP

HAXBY ROAD • YORK • ENGLAND

Telegrams: Coordinate, York

Telephone: York 24112

LONDON ADDRESS:

Vickers House, Millbank, London, S.W.1.

Telegrams: Coordinate, London, S.W.1.

Telephone: TATe Gallery 7777

**Technical sales representatives operate throughout
Great Britain**

**OFFICES
ABROAD**

CANADA

TORONTO 23 Rainside Road, Don Mills, Ontario

Telephone: Hlckory 7-5135 Telegrams: Coordinate, Toronto

MONTREAL 295 Villeneuve Street West, Montreal, Quebec

Telephone: Crescent 6-4863

EDMONTON 10454 82nd Avenue, Edmonton, Alberta

Telephone: GENEva 9-3324

OTTAWA 1396 Merivale Road, Ottawa, Ontario

Telephone: Parkway 83881

S. RHODESIA

SALISBURY 710-711 Linqenda House, Baker Avenue, Salisbury, C.I.

P.O. Box 3777. Telephone: 28814. Telegrams: Coordinate, Salisbury, Rhodesia

S. AFRICA

COOKE, TROUGHTON & SIMMS SOUTH AFRICA (PTY) LTD.

JOHANNESBURG M.B.S. Building, 16 Wolmarans Street, P.O. Box 7131

Telephone: 22-1394 Telegrams: Coordinate, Johannesburg

CAPE TOWN Sun Building, St. Georges Street, P.O. Box 1552

Telephone: 41-2135 Telegrams: Coordinate, Cape Town

U.S.A.

COOKE, TROUGHTON & SIMMS INC.

MALDEN 91 Waite Street, Malden 48, Massachusetts, P.O. Box 93 (nr. Boston)

Telephone: DAVenport 4-6666 Telegrams: Coordinate, Boston



COOKE TROUGHTON & SIMMS LTD.
YORK • ENGLAND

VICKERS INSTRUMENTS LTD.

Successors to: Cooke Troughton & Simms Ltd.—C. Baker Instruments Ltd.

HAXBY ROAD, YORK

Telephone: York 24112

Telegrams: Coordinate York

PURLEY WAY, CROYDON

Telephone: Croydon 3845

Telegrams: Optivorum Croydon

VICKERS M55 MICROSCOPE

JULY 1964

		Price			Number Required	Total		
		£	s.	d.		£	s.	d.
Vickers Fifty-Five Microscope, comprising:— microscope block, main casing, incident illuminator unit, lamp slide-way, lamp condenser unit, zoom projection unit, small projection mirror, large pro- jection mirror, xenon lamp, electrical equipment, end panel normal illumination, 2 in. square, ground glass diffuser; 2 in. square, heat absorbing filter (ON.20); water trough, plastic instrument cover, stop collar for 10× Kellner eyepiece, 10× Kellner eyepiece with graticule, objective hole plugs, coarse focusing adjusting spanner, instruction book, goggles for xenon lamp, frame for focusing screen, and fixed wiring for integrator	M550000	<input type="checkbox"/>	1140	0	0			
Gliding stage with joystick control	M551515	<input type="checkbox"/>	44	0	0			
Revolving sextuple objective changer	M550560	<input type="checkbox"/>	19	4	0			
Magnification changer unit	M550011	<input type="checkbox"/>	100	16	0	X		
Monocular head	M550670	<input type="checkbox"/>	4	4	0			
Plain front panel	M551955	<input type="checkbox"/>	2	8	0			
Double plate holder for ½-plates	M505647	<input type="checkbox"/>	7	12	0			
Two ¼-plate adaptors for M505647	M505960	<input type="checkbox"/>	4	16	0			
Micro crystalline wax focusing screen	M551700	<input type="checkbox"/>	4	16	0			
Accessory box for objectives and eyepieces	M552310	<input type="checkbox"/>	3	18	0			
M55 Microscope, as above, for manually operated photography	M550001	<input type="checkbox"/>	1331	14	0			

		Price			Number Required	Total		
		£	s.	d.		£	s.	d.
Vickers Fifty-Five Microscope, comprising:— microscope block, main casing, incident illuminator unit, lamp slide-way, lamp condenser unit, zoom projection unit, small projection mirror, large pro- jection mirror, xenon lamp, electrical equipment, end panel normal illumination, 2 in. square, ground glass diffuser; 2 in. square, heat absorbing filter (ON.20); water trough, plastic instrument cover, stop collar for 10× Kellner eyepiece, 10× Kellner eyepiece with graticule, objective hole plugs, coarse focusing adjusting spanner, instruction book, goggles for xenon lamp, frame for focusing screen, and fixed wiring for integrator	M550000	<input type="checkbox"/>	1140	0	0			
Swing-out transmitted light bracket	M550500	<input type="checkbox"/>	88	0	0			
Centring, rotating and gliding stage	M550300	<input type="checkbox"/>	120	0	0			
Revolving sextuple objective changer	M550560	<input type="checkbox"/>	19	4	0			
Revolving quintuple objective changer	M551545	<input type="checkbox"/>	24	0	0			
Magnification changer unit with analyser unit	M550012	<input type="checkbox"/>	129	12	0	<input checked="" type="checkbox"/>		
Binocular head	M550680	<input type="checkbox"/>	40	0	0			
Motorised 35 mm. camera unit	M550015	<input type="checkbox"/>	160	0	0	<input checked="" type="checkbox"/>		
Automatic integrating photographic timing unit	M550017	<input type="checkbox"/>	282	0	0			
Double plate holder for ½-plates	M505647	<input type="checkbox"/>	7	12	0			
Two ¼-plate adaptors for M505647	M505960	<input type="checkbox"/>	4	16	0			
Micro crystalline wax focusing screen	M551700	<input type="checkbox"/>	4	16	0			
Desk for microscope	M552380	<input type="checkbox"/>	98	0	0			
M55 Microscope, as above, for automatic 35 mm. and plate photography		M550002	<input type="checkbox"/>	2118	0	0		2118 0 0

PHOTOGRAPHIC EQUIPMENT

Automatic integrating photographic timer (including motorised focal plane shutter M551850)	M550017	<input type="checkbox"/>	282	0	0			
Motorised focal plane shutter	M551850	<input type="checkbox"/>	122	0	0			
Motorised automatic 35 mm. camera	M550015	<input type="checkbox"/>	160	0	0			
Manual control 35 mm. camera unit	M550014	<input type="checkbox"/>	97	18	0			
Additional laboratory cassette	M551854	<input type="checkbox"/>	67	4	0			
MPP camera back (to take 5 in. × 4 in. plates), wood adapter, focusing screen and hood	M552095	<input type="checkbox"/>	19	16	0			
MPP camera back, wood adapter, focusing screen, hood and polaroid land cut film holder (5 in. × 4 in.)	M552100	<input type="checkbox"/>	60	8	0			
MPP camera back, wood adapter, focusing hood and polaroid J.66 roll film back	M552105	<input type="checkbox"/>	90	2	0			
Simple photographic shutter	M505705	<input type="checkbox"/>	11	14	0			

		Price	Number		Total	
		£ s. d.	Required	£ s. d.	£ s. d.	
Screen focusing magnifier	M505890 <input type="checkbox"/>	8 0 0				
Double plate holder for 7 in. × 5 in. plates	M505664 <input type="checkbox"/>	9 12 0				
Two adapters for 5 in. × 4 in. plate for M505664	M505689 <input type="checkbox"/>	3 4 0				
Double plate holder for half plates (6½ in. × 4¾ in.)	M505647 <input type="checkbox"/>	7 12 0				
Two adapters for ¼-plates (3¼ in. × 4¼ in.) for M505647	M505960 <input type="checkbox"/>	4 16 0				
Two adapters for 9 cms. × 12 cms. plates for M505647	M505790 <input type="checkbox"/>	3 16 0				
7 in. × 5 in. cut film holder for M505664	M505931 <input type="checkbox"/>	1 14 0				
5 in. × 4 in. cut film holder for M505689	M505932 <input type="checkbox"/>	16 0				
6½ in. × 4¾ in. cut film holder for M505647	M505933 <input type="checkbox"/>	10 0				
3¼ in. × 4¼ in. cut film holder for M.505960	M505934 <input type="checkbox"/>	10 0				
Frame for focusing screens	M551610 <input type="checkbox"/>	2 16 0				
Focusing screen ruled with cross lines	M505904 <input type="checkbox"/>	3 12 0				
Focusing screen with horizontal and vertical lines intersecting at centre of the plate, each graduated every 100 mms	M505906 <input type="checkbox"/>	3 10 0				
Focusing screen with grid 10 cms. × 10 cms. ruled in mms. with break at intersections of 1 cm. lines	M505907 <input type="checkbox"/>	3 14 0				
Focusing screen with grid 10 cms. × 10 cms. ruled in cms.	M505908 <input type="checkbox"/>	3 18 0				
STAGES						
Divided rotating stage with micrometer movements and gliding top plate	M550300 <input type="checkbox"/>	120 0 0				
Gliding stage with joystick control	M551515 <input type="checkbox"/>	44 0 0				
VIEWING HEADS						
Binocular head	M550680 <input type="checkbox"/>	40 0 0				
Monocular head	M550670 <input type="checkbox"/>	4 4 0				
MAGNIFICATION CHANGERS						
Magnification changer with focusing Bertrand lens	M550011 <input type="checkbox"/>	100 16 0				
Magnification changer with focusing Bertrand lens and rotating polarizer, rotating quartz plate (both swing out) and rotating compensator slot	M550012 <input type="checkbox"/>	129 12 0				
OBJECTIVE CHANGERS						
Revolving sextuple objective changer for incident light field objectives	M550560 <input type="checkbox"/>	19 4 0				
*Revolving triple objective changer for incident dark ground objectives with catoptric condensers and for use with micro-hardness testing objectives	M552145 <input type="checkbox"/>	16 8 0				
<i>*The quintuple objective changer, M550550, previously supplied is now obsolete</i>						

		Price			Number Required	Total		
		£	s.	d.		£	s.	d.
Revolving quintuple objective changer incorporating a 2× correcting lens for use with standard 160 mm. tube length, transmitted light objectives	M551545 <input type="checkbox"/>	24	0	0				

*INCIDENT LIGHT PHASE CONTRAST

Incident phase contrast unit for positive and negative phase contrast	M550570 <input type="checkbox"/>	176	0	0			
---	----------------------------------	-----	---	---	--	--	--

* This apparatus uses the standard ∞ corrected objectives

TRANSMITTED LIGHT PHASE CONTRAST

Swing out substage bracket with focusing substage and mirror	M550500 <input type="checkbox"/>	88	0	0			
--	----------------------------------	----	---	---	--	--	--

Phase contrast condenser with rotating disc of four substage annuli	M555625 <input type="checkbox"/>	56	0	0			
---	----------------------------------	----	---	---	--	--	--

Long working distance phase contrast condenser with rotating disc of four substage annuli	M555642 <input type="checkbox"/>	52	0	0			
---	----------------------------------	----	---	---	--	--	--

(Phase contrast objectives are listed under objectives)

ILLUMINATION

High power tungsten filament lamp in lamp house with condenser (for transmitted light only)	M550013 <input type="checkbox"/>	62	16	0			
---	----------------------------------	----	----	---	--	--	--

Spare xenon bulb CSX 150 W.	M550828 <input type="checkbox"/>	21	0	0			
-----------------------------	----------------------------------	----	---	---	--	--	--

Spare high pressure mercury bulb	M551649 <input type="checkbox"/>	18	10	0			
----------------------------------	----------------------------------	----	----	---	--	--	--

Spare 6V., 30W. filament bulb	E855 <input type="checkbox"/>		14	8			
-------------------------------	-------------------------------	--	----	---	--	--	--

62 16

MICRO HARDNESS TESTING

Micro hardness testing equipment, including transmitter, indenter objective in centring mount, 50× measuring objective in centring mount, revolving triple objective changer, monocular head with centring filar micrometer eyepiece, specimen holder and wood box	M550016 <input type="checkbox"/>	200	0	0			
--	----------------------------------	-----	---	---	--	--	--

OBLIQUE ILLUMINATION UNIT

Oblique illumination unit with lens to go into incident illuminator tube	M550020 <input type="checkbox"/>	26	0	0			
--	----------------------------------	----	---	---	--	--	--

MACRO EQUIPMENT

Macro Objectives				Price		
Focal length		Magnification		£	s.	d.
5.0 in.	f/4.5	5×	M551325 <input type="checkbox"/>	15	8	0
2.6 in.	f/4.5	10×	M551340 <input type="checkbox"/>	8	16	0
1.9 in.	f/3.5	15×	M551360 <input type="checkbox"/>	12	16	0

1250

For Macro Incident Illumination

		£	Price s.	d.		
*Macro base unit with revolving changer for macro objectives	M551270	15	10	0	1	
Incident illuminator for 10× and 15× macro objectives	M551795	19	4	0	1	
Incident illuminator for 5× macro objective	M551815	16	16	0	1	
Projection lens for incident light macro with 5× objective	M551830	4	12	0	1	
Mirror for oblique illuminator	M551310	3	12	0	1	

For Macro Transmitted Illumination

*Macro base unit with revolving changer for macro objectives	M551270	15	10	0		
Swing-out mirror bracket	M550500	88	0	0	1	
Auxiliary condenser for transmitted light macro	M551825	5	12	0	1	
Projection lens for transmitted light macro with 5×, 10× and 15× objectives	M551990	4	4	0	1	

**The macro base unit M551270 is common to both incident and transmitted illumination*

POLARIZING ACCESSORIES

Monocular head with swing out focusing Bertrand lens	M551897	48	0	0	1	
Polarizing substage attachment	M551040	24	8	0	1	
Graduated polarizing cap (for incident light only)	M552065	6	8	0	1	

(Strain free objectives are listed separately)

Compensators

Quartz wedge (not graduated)	M552075	7	2	0	1	
Quartz wedge (graduated)	M552080	9	16	0	1	
Mica ¼-wave plate	M552085	2	8	0	1	
de Sénarmont compensator	M552090	3	4	0	1	
Nakamura half shadow plate	M552056	9	12	0	1	
Elliptic compensator	M550950	9	12	0	1	

ACCESSORIES

Small bottle of non-drying immersion oil as supplied with each immersion objective (ALP ₁) nd. _D =1.524	M001195	4	0	0		
8 oz. bottle of above oil	M001196	2	0	0		
Small bottle of non-drying immersion oil (ALP ₂) nd. _D =1.515	M001197	4	0	0		
8 oz. bottle of above oil	M001198	2	0	0		
Eyepiece micrometer 1 cm. divided into 100 parts	M001581	2	12	0		

			Price			Number Required	Total		
			£	s.	d.		£	s.	d.
Glass stage micrometer 1 mm. divided into 100 parts	M001586	<input type="checkbox"/>	4	0	0				
Polished metal stage micrometer 1 mm. divided into 100 parts	M001591	<input type="checkbox"/>	5	12	0				
Filar micrometer eyepiece 10× normal field, reading directly to 0.01 mm. on a micrometer drum with estimation to 0.002 mm.	M001555	<input type="checkbox"/>	32	8	0				
As above with 15× eyepiece	M001556	<input type="checkbox"/>	32	8	0				
Water trough	M505970	<input type="checkbox"/>	1	4	0				
Heat absorbing filter (O.N. 20)	M505608	<input type="checkbox"/>	1	0	0				
Soft plastic dust cover (supplied with the instrument)	M552300	<input type="checkbox"/>	1	14	0				

EYEPIECES

Huyghens

	Old No.	New No.						
6× single	M001506	M040100	<input type="checkbox"/>	1	14	0		
8× single	M001511	M040700	<input type="checkbox"/>	1	14	0		
10× single	M001516	M040300	<input type="checkbox"/>	1	14	0		
6× pair	M001507	M040120	<input type="checkbox"/>	3	14	0		
8× pair	M001512	M040720	<input type="checkbox"/>	3	14	0		
10× pair	M001517	M040320	<input type="checkbox"/>	3	14	0		

Compensating

6× single	M001526	M041100	<input type="checkbox"/>	2	12	0		
8× single	M001531	M041700	<input type="checkbox"/>	2	16	0		
10× single		M041300	<input type="checkbox"/>	3	8	0		
20× single	M001551	M041602	<input type="checkbox"/>	5	16	0		
6× pair	M001527	M041120	<input type="checkbox"/>	5	10	0		
8× pair	M001532	M041720	<input type="checkbox"/>	5	18	0		
10× pair		M041320	<input type="checkbox"/>	7	2	0		
20× pair	M001552	M041622	<input type="checkbox"/>	11	18	0		

Kellner

10× single	M001541	M042302	<input type="checkbox"/>	4	14	0		
10× pair	M001542	M042322	<input type="checkbox"/>	9	14	0		

Complan

10× single	M019536	M041301	<input type="checkbox"/>	7	2	0		
10× pair	M019537	M041321	<input type="checkbox"/>	14	10	0		

				Price			Number Required	Total		
				£	s.	d.		£	s.	d.
CONDENSERS FOR TRANSMITTED LIGHT										
Centring condenser mount			M001376	<input type="checkbox"/>	7	0	0			
Abbe condenser			M001382	<input type="checkbox"/>	3	8	0			
Aplanatic condenser			M001383	<input type="checkbox"/>	7	16	0	X		
Achromatic condenser (N.A. 1.30 when oil immersed)			M001391	<input type="checkbox"/>	14	16	0			
Focusing dark ground condenser			M001396	<input type="checkbox"/>	22	8	0	X		
Low power dark ground condenser with fixed focus N.A. 0.7-0.8			M001394	<input type="checkbox"/>	18	0	0			
STRAIN FREE CONDENSERS FOR TRANSMITTED LIGHT										
Aplanatic condenser			M007884	<input type="checkbox"/>	7	16	0			
Achromatic condenser (N.A. 1.30 when oil immersed)			M007891	<input type="checkbox"/>	14	16	0			
OBJECTIVES FOR UNCOVERED SPECIMENS										
Infinite tube length										
Achromatic										
	<i>Magnification</i>	<i>N.A.</i>								
	3.5×	0.05	M023052	<input type="checkbox"/>	2	16	0			
	6×	0.10	M022052	<input type="checkbox"/>	2	18	0			
	10×	0.15	M022152	<input type="checkbox"/>	5	12	0			
	15×	0.25	M022352	<input type="checkbox"/>	5	16	0			
	30×	0.5	M022452	<input type="checkbox"/>	8	18	0			
	50×	0.8	M022952	<input type="checkbox"/>	9	0	0			
	85×	0.85	M023152	<input type="checkbox"/>	11	6	0			
	140× oil	1.30	M022652	<input type="checkbox"/>	14	0	0			
Fluorite										
	140× oil	1.30	M023552	<input type="checkbox"/>	32	0	0			
Microplan										
	30×	0.5	M025452	<input type="checkbox"/>	25	0	0			
	40×	0.65	M025352	<input type="checkbox"/>	34	4	0			
Apochromatic										
	17×	0.3	M024052	<input type="checkbox"/>	27	0	0			
	50×	0.95	M023852	<input type="checkbox"/>	53	0	0			
	85×	0.95	M023352	<input type="checkbox"/>	44	12	0			
	115× oil	1.32	M023752	<input type="checkbox"/>	43	0	0			
Funnel stop for objectives with N.A. above 1.0 when used with M001396 dark ground condenser			M001397	<input type="checkbox"/>	6	0	0			

OBJECTIVES FOR UNCOVERED SPECIMENS (continued)

Infinite tube length

Strain Free for Polarized Light

Achromatic

Magnification N.A.

Magnification	N.A.		Price £ s. d.
3.5 ×	0.05	M023054 <input type="checkbox"/>	3 6 0
6 ×	0.10	M022054 <input type="checkbox"/>	3 6 0
10 ×	0.15	M022154 <input type="checkbox"/>	6 0 0
15 ×	0.25	M022354 <input type="checkbox"/>	5 16 0
30 ×	0.5	M022454 <input type="checkbox"/>	9 4 0
50 ×	0.8	M022954 <input type="checkbox"/>	9 10 0
85 ×	0.85	M023154 <input type="checkbox"/>	11 8 0
140 × oil	1.3	M022654 <input type="checkbox"/>	14 2 0

For Incident Light, Dark Field

15 ×	0.25	M023252 <input type="checkbox"/>	10 12 0
30 ×	0.5	M023452 <input type="checkbox"/>	9 8 0
50 ×	0.65	M024152 <input type="checkbox"/>	10 2 0

Catoptric Condensers for Incident Light Dark Field

Catoptric condenser for objectives M023452 or M024152	M551896 <input type="checkbox"/>	8 16 0
Catoptric condenser for objective M023252	M551895 <input type="checkbox"/>	6 16 0

OBJECTIVES FOR COVERED SPECIMENS

160 mm. Tube Length to be used with objective changer M551545 having a 2 × corrector lens

Achromatic

3 ×	0.10	M022011 <input type="checkbox"/>	2 14 0
5 ×	0.15	M022101 <input type="checkbox"/>	5 18 0
10 ×	0.25	M022301 <input type="checkbox"/>	4 10 0
20 ×	0.50	M022401 <input type="checkbox"/>	8 16 0
40 ×	0.65	M022501 <input type="checkbox"/>	7 8 0
40 ×	0.85	M022901 <input type="checkbox"/>	10 10 0
95 × oil	1.30	M022601 <input type="checkbox"/>	11 6 0

Fluorite

95 × oil	1.30	M023501 <input type="checkbox"/>	31 10 0
----------	------	----------------------------------	---------

Microplan

10 ×	0.25	M025111 <input type="checkbox"/>	9 16 0
20 ×	0.50	M025411 <input type="checkbox"/>	15 14 0
40 ×	0.70	M025211 <input type="checkbox"/>	12 18 0

Apochromatic

10 ×	0.30	M024001 <input type="checkbox"/>	27 0 0
40 ×	0.95	M023801 <input type="checkbox"/>	53 0 0
80 ×	1.32	M023701 <input type="checkbox"/>	43 0 0

OBJECTIVES FOR COVERED SPECIMENS (continued)

Objectives, Strain Free for Polarized Light

Achromatic

			Price			Number Required	Total		
			£	s.	d.		£	s.	d.
3 ×	0.10	M022013	3	4	0				
5 ×	0.15	M022103	6	6	0				
10 ×	0.25	M022303	4	18	0				
20 ×	0.50	M022403	9	2	0				
40 ×	0.65	M022503	7	12	0				
40 ×	0.85	M022903	10	12	0				
95 × oil	1.30	M022603	11	8	0				

Objectives for Transmitted Light Phase Contrast

Achromatic

10 ×	0.25	M022205	5	18	0			
20 ×	0.5	M022405	9	14	0			
40 ×	0.65	M022505	9	4	0			
95 × oil	1.30	M022605	13	0	0			
Fluorite								
45 × oil	0.95	M023605	36	10	0			

DELIVERY

Packing and carriage

TOTAL

Yours faithfully,
for VICKERS INSTRUMENTS LTD.

Signed

SUMMARY in numerical order

	£	s.	d.		£	s.	d.		£	s.	d.
E855		14	8	M022652	14	0	0	M505960	4	16	0
M001195		4	0	M022654	14	2	0	M505970	1	4	0
M001196	2	0	0	M022901	10	10	0	M550000	1140	0	0
M001197		4	0	M022903	10	12	0	M550001	1331	14	0
M001198	2	0	0	M022952	9	0	0	M550002	2118	0	0
M001376	7	0	0	M022954	9	10	0	M550011	100	16	0
M001382	3	8	0	M023052	2	16	0	M550012	129	12	0
M001383	7	16	0	M023054	3	6	0	M550013	62	16	0
M001391	14	16	0	M023152	11	6	0	M550014	97	18	0
M001394	18	0	0	M023154	11	8	0	M550015	160	0	0
M001396	22	8	0	M023252	10	12	0	M550016	200	0	0
M001397		6	0	M023352	44	12	0	M550017	282	0	0
M001506	1	14	0	M023452	9	8	0	M550020	26	0	0
M001507	3	14	0	M023501	31	10	0	M550300	120	0	0
M001511	1	14	0	M023552	32	0	0	M550500	88	0	0
M001512	3	14	0	M023605	36	10	0	M550560	19	4	0
M001516	1	14	0	M023701	43	0	0	M550570	176	0	0
M001517	3	14	0	M023752	43	0	0	M550670	4	4	0
M001526	2	12	0	M023801	53	0	0	M550680	40	0	0
M001527	5	10	0	M023852	53	0	0	M550828	21	0	0
M001531	2	16	0	M024001	27	0	0	M550950	9	12	0
M001532	5	18	0	M024052	27	0	0	M551040	24	8	0
M001541	4	14	0	M024152	10	2	0	M551270	15	10	0
M001542	9	14	0	M025111	9	16	0	M551310	3	12	0
M001551	5	16	0	M025211	12	18	0	M551325	15	8	0
M001552	11	18	0	M025352	34	4	0	M551340	8	16	0
M001555	32	8	0	M025411	15	14	0	M551360	12	16	0
M001556	32	8	0	M025452	25	0	0	M551515	44	0	0
M001581	2	12	0	M040100	1	14	0	M551545	24	0	0
M001586	4	0	0	M040120	3	14	0	M551610	2	16	0
M001591	5	12	0	M040300	1	14	0	M551649	18	10	0
M007884	7	16	0	M040320	3	14	0	M551700	4	16	0
M007891	14	16	0	M040700	1	14	0	M551795	19	4	0
M019536	7	2	0	M040720	3	14	0	M551815	16	16	0
M019537	14	10	0	M041100	2	12	0	M551825	5	12	0
M022011	2	14	0	M041120	5	10	0	M551830	4	12	0
M022013	3	4	0	M041300	3	8	0	M551850	122	0	0
M022052	2	18	0	M041320	7	2	0	M551854	67	4	0
M022054	3	6	0	M041602	5	16	0	M551895	6	16	0
M022101	5	18	0	M041622	11	18	0	M551896	8	16	0
M022103	6	6	0	M041700	2	16	0	M551897	48	0	0
M022152	5	12	0	M041720	5	18	0	M551955	2	8	0
M022154	6	0	0	M042302	4	14	0	M551990	4	4	0
M022205	5	18	0	M042322	9	14	0	M552056	9	12	0
M022301	4	10	0	M505608	1	0	0	M552065	6	8	0
M022303	4	18	0	M505647	7	12	0	M552075	7	2	0
M022352	5	16	0	M505664	9	12	0	M552080	9	16	0
M022354	5	16	0	M505689	3	4	0	M552085	2	8	0
M022401	8	16	0	M505705	11	14	0	M552090	3	4	0
M022403	9	2	0	M505790	3	16	0	M552095	19	16	0
M022405	9	14	0	M505890	8	0	0	M552100	60	8	0
M022452	8	18	0	M505904	3	12	0	M552105	90	2	0
M022454	9	4	0	M505906	3	10	0	M552145	16	8	0
M022501	7	8	0	M505907	3	14	0	M552300	1	14	0
M022503	7	12	0	M505908	3	18	0	M552310	3	18	0
M022505	9	4	0	M505931	1	14	0	M552380	98	0	0
M022601	11	6	0	M505932		16	0	M555625	56	0	0
M022603	11	8	0	M505933		10	0	M555642	52	0	0
M022605	13	0	0	M505934		10	0				